

What's New

April 2008

GERBER CHAIR MATES NEWS

New Item

Introducing

THE SIDEKICK

Many people enjoy target shooting and hunting.

The same is true for people confined to a wheelchair. The enjoyment and pride of such hobbies can be particularly important to a person confined to a wheel chair as it provides a challenging and interesting activity in addition to motivating the individual to venture outdoors and become active.

Recognizing this need, Gerber Chair Mates developed The Sidekick in order to provide the BEST in comfort and support for those individuals that hunt or target shoot from a wheelchair.

The Sidekick mounts easily onto the arms of the wheelchair with hook & loop strips, which allows for ease in release for entry and exit of the chair.

The pad can be adjusted to the right or left for target practice and is a comfortable place for the hunter to rest their gun, crossbow or arm for a more steady shot.

The support pad can be sized according to individual need and it adjusts vertically for height. The standard framing is available in Black or Olive Drab. Camouflage can be done at an additional cost.

One of the benefits of using The Sidekick is how it aids in facilitating proper posture when shooting to minimize injury. Additionally, it allows the individual to continue to participate in outdoor hobbies they enjoy while promoting a feeling of independence.

The Sidekick is compatible with most wheelchairs on the market.

Visit our New Products page for ordering information.



Ki Mobility Wheelchairs

NEW NEW NEW NEW NEW NEW NEW NEW NEW NEW NEW

Over the past several months, a few customers have inquired about the Gerber Swing-Away Stump Support fitting the Ki Mobility Wheelchair.

Due to a difference in the front rigging on these chairs, our current models were not adaptable.

In order to stay current with these chairs, we contacted the company and are currently working with them to develop an adapter that will allow our support to attach to these chairs.

We will keep you updated on our progress.



Rehabilitation Facilities

Q. Why should you select a Gerber Swing-Away Stump Support for your amputee patients?

A. The Gerber Swing-Away Stump Support is specifically designed to hold the residual limb in the proper healing position, and can be attached to the wheelchair without tools.

Doctors and Therapists needing to care for the patient while in the chair can easily swing it out of the way. It doesn't require a service call to a technician for removal when/if the patient wants the leg rests on their chair. It's easy enough that even the patient can do it.

MANUFACTURED HERE

There are several amputee stump supports on the market today. The Gerber stump support offers an easy to assemble unit with Swing-Away action. It is small and lightweight. Our stump support is made in the USA so there is minimal chance of factory recalls and inferior manufacturing.



Let us help

If you experience any problem with the Swing-Away Stump Support adapting to a particular chair, please contact us immediately. Wheelchair manufacturers change styles, models and designs, often without our knowledge. By bringing it to our attention, we can work together to solve the problem. Email us at sales@gerberchairmates.com