

Gerber Chair Mates, Inc

Newsletter

October 2008

KI Mobility Compatibility

After several months of work, Gerber Chair Mates has designed a Swing-Away Stump Support to fit the KI Mobility Wheelchairs. The design has the same features and benefits as the other Gerber Models

The model numbers for the supports are GCM-50IR for right side and GCM-502L for left side. The retail cost for the support is \$340.00. It maintains the same HCPCS Code EI020.

Gerber Chair Mates is grateful to KI Mobility owner, Doug Munsey for the assistance he provided. His willingness to offer input, measurements and a frame gave our manufacturing team the basis needed to develop the support.

The supports are in stock and ready for shipment.

Adjustable Support

Several customers have inquired whether the Gerber Stump Support is adjustable on an angle, vertically. Typically our manufacturing department would custom make the "L" tube at the angle requested by the therapist. However, this does leave it at a fixed angle. Realizing that patients may need to adjust the support at different angles for smoother transitions, our manufacturing dept has begun work on having the Gerber Support adjust to various angles.

We will keep you informed on our progress.

Manufactured in the USA

There are several amputee stump supports on the market today. The Gerber Stump Support offers an easy to assemble unit with Swing-Away action. It is small and lightweight. Our stump support is made in the [USA](#) so there is minimal chance of factory recalls and inferior manufacturing.



Let us help

If you experience any problem with the Swing-Away Stump Support adapting to a particular chair, please contact us immediately. Wheelchair manufacturers change styles, models and designs, often without our knowledge. By bringing it to our attention, we can work together to solve the problem. Email us at sales@gerberchairmates.com